



Performance And Value

U.S. Manufactured, Ceramic Bearing, MULTIFlex* Compatible Highspeed Handpieces

- Ceramic bearings for extended life and smooth operation
- Low maintenance required

Choose The Quest Head Size To Fit Your Needs



100 Series

Exceptional Visibility and Accessibility

- Compact head for exceptional visibility and accessibility
- Fast cutting
- Easy posterior access
- Twin beam optics (optional)
- Secure push button bur changing
- Chairside turbine and chuck replacement
- Autoclavable and chemiclavable
- Ergonomic, hygienic grip

Heads shown
actual size



300 Series

Smooth and Powerful

- Powerful torque
- Heavy duty construction
- Smooth, quiet operation
- Twin beam optics (optional)
- Secure push button bur changing
- Chairside turbine and chuck replacement
- Autoclavable and chemiclavable
- Ergonomic, hygienic grip



Handpieces & Coupler shown
actual size

Quest 150

Power and durability
in a compact head
with stainless steel
construction.
50 grams



Quest 350

Triple port water spray
and stainless steel
construction with power
for your toughest
cutting requirements.
56 grams



Quest is also
available in
standard four
and five-hole
models.



*MULTIFlex is a trademark
of KaVo

Quest

handpieces

All Quest Handpieces are sold exclusively through your dental dealer. Quest handpieces are manufactured in the U.S.A. by Lares Research with a tradition of 55 years of manufacturing excellence. This combination provides you, the Quest owner, with solid product manufacturing experience matched with the local service of your dealer.

Quest models use stainless steel construction for easy to clean, hygienic finishes and low maintenance ceramic bearings for long life. Quest handpieces offer several optional features; two different head sizes, fiber optics and swivel capability that is compatible with the popular MULTIflex* coupler.

Connection

Quest models offer the MULTIflex* alternative you can afford



Optics

Optional twin beam fiber optics provide enhanced shadow-free illumination

Cooling

Quest 300 models feature multi-port water spray



Specifications

	100 Series	300 Series	Weight
Speed	450,000 rpm	360,000 rpm	Quest 130/150 50g
Power	15.3 watts	20.8 watts	Quest 330/350 56g
Minimum air pressure	32 psi/220 kPa	32 psi/220 kPa	Quest 110/120 41g
Maximum air pressure	40 psi/275 kPa	40 psi/275 kPa	Quest 310/330 48g
Average overall head length (with 19 mm bur inserted)	21.5 mm	22 mm	
Head length	11.9 mm	14.65 mm	
Head diameter	10.2 mm	13 mm	

Quest handpieces provide you the tool to meet your practice needs, delivered by your favorite dealer.

Your Local Quest Dealer:



Made in the U.S.A.

295 Lockheed Avenue
Chico, CA 95973